



SCHILLER
The Art of Diagnostics



CARDIOVIT AT-102 G2

Less
6%



R 86 150.00 Excl VAT - Unit only

Less
8%



R 101 400.00 Excl VAT —

Unit Including Mobile trolley

Easy, smart and state-of-the-art

Best signal quality and diagnostic algorithms for adult and paediatric patients

Standard Wi-Fi and LAN connectivity, combined with ECG review on the display and a large memory, supports paperless workflows and cost saving

Schiller Link application to save recordings in PDF format on PC

Spirometry and basic exercise ECG optionally available

Sealed alphanumeric keyboard, easy to clean and highly hygienic

Navigation keys for high efficiency in any workflow

Internal memory: up to 350 ECG, 100 Resting Rhythm Recordings and 10 exercise ECG

Available units : 15

Price valid while stock last



SCHILLER

The Art of Diagnostics



CARDIOVIT MS-2015

Less
38%



R 66 250.00 Excl VAT Including
Mobile trolley

- Advanced 12/16-channel ECG system with touch screen technology.
- Resting ECG, manual printout, resting rhythm up to 4 minutes
- Colour 15" touch screen
- Battery operation
- Large internal memory
- On screen review of ECG
- High resolution thermal printer (paper size A4)
- ECG measurements
- PDF export on USB
- Interfaces: LAN, USB
- 16-Lead ECG optional with 14-Lead ECG cable

Available units : 3

Price valid while stock last



SCHILLER
The Art of Diagnostics



CARDIOVIT FT-1



R 46 900.00 Excl VAT

- Maximum performance in a compact electrocardiograph
- For users who value state-of-the art technology
- Power and flexibility in a portable ECG - weighing just over 1 Kg
- Fast and secure bidirectional Wi-Fi communication, SCHILLER LINK: Easy-to-use software to import patient data from the EMR system and export customized PDF reports directly on the PC.
- Outstanding signal quality for adult and paediatric ECG
- 8" Inch high resolution multi-touch screen
- PAEDIATRIC ECG: Thanks to very high sampling frequency and a large frequency band, offers optimal signal quality which makes it the ideal tool for paediatric ECG.
- Anatomical hook up adviser helps to place the electrodes correctly with zoom-in and rotate functions for optimal view

Available units : 1

Price valid while stock last



SCHILLER

The Art of Diagnostics



***BR-102 plus* Ambulatory Blood Pressure**

Less
24%



R 38 250.00 Excl VAT

- BR-102 plus Blood pressure recorder
- The new standard in ambulatory blood pressure monitoring.
- The Third Generation of BR-102 plus sets the standard in ABPM by meeting the highest clinical standards, optimizing patient comfort through innovation and giving clinicians practical
- interpretive analysis and reporting tools.
- Approved by BHS : ESH : AAMI: dabl
- Standard unit including:
- Auscultatory measurement method (main) in combination with oscillometric backup measurement (Riva-Rocci, Korotkoff / Oscillometric back-up measurement)
- Darwin2 ABPM evaluation software
- Automated analysis

Available units : 1

Price valid while stock last



SCHILLER

The Art of Diagnostics



***BR-102 plus PWA* Ambulatory Blood Pressure with Pulse wave analysis**

Less
48%



R 43 600.00 Excl VAT

- Pulse Wave Analysis (PWA), a new method for the early detection of cardiovascular diseases, is increasingly attracting attention.
- Pulse wave analysis can be performed in a 24-hours cycle to receive a full day PWA profile including changes during the night period or, for increased patient comfort, a spot measurement of PWA can be performed only at the very beginning of each 24-hours cycle.
- Using the new spot measurement option, one can analyse pulse wave also in only a 20 minutes cycle.
- Standard unit including:
 - Darwin2 ABPM evaluation software
 - Auscultatory measurement method (main) in combination with oscillometric backup measurement (Riva-Rocci, Korotkoff / Oscillometric back-up measurement)
 - Approved by BHS : ESH : AAMI: dabl

Available units : 1

Price valid while stock last



SCHILLER

The Art of Diagnostics



Spiroscout with LFX



R 40 400.00 Excl VAT including
3L Calibration syringe and
Weather station. Computer
not included

- SpiroScout incorporates an ultrasonic flow sensor to measure the flow of air in and out of the patients lungs. Ultrasonic flow measurement eliminates problems associated with traditional methods of flow measurement. There are no moving parts, no filters to change, no screens to catch sputum, and no disposables to calibrate. It is calibration free.
- The new LFX PFT software is the highest quality advanced spirometry software you can buy. It is simple to use and employs a state of the art multi-window, graphical interface. Up to 41 Spirometry indices can be measured, previous test curves overlaid for comparison and both open and closed Flow/Volume Loop test techniques used. It also features an animated child incentive device and the long term trending of patient results.
- Computer or laptop required for operation, Windows 10

Available units : 4

Price valid while stock last



SCHILLER

The Art of Diagnostics



FRED EASY AED

Less
33%



R 26 485.00 Excl VAT

- FRED easy semi-automatic, public access defibrillator with Biphasic Multipulse Biowave technology
- Adapted voice and displayed instructions for first aid rescuers together with a fast rhythm analysis and shock delivery guide the rescuer during the resuscitation process and minimise the hands-off time.
- Works with disposable LiMnO2 cell

Available units : 3

Price valid while stock last



SCHILLER

The Art of Diagnostics



FRED EASY AED TRAINER

Less
42%



- FRED EASY TRAINER
- Allows training in Semi-automatic or in Fully-Automatic mode

R 5 890.00 Excl VAT

Available units : 3

Price valid while stock last



SCHILLER

The Art of Diagnostics



Fred Easyport AED

Less
55%



- Semi-automatic pocket defibrillator with Biphasic Multipulse® Biowave technology
- The worlds first pocket defibrillator
- High resolution LCD
- Easy to operate
- Quick charging time
- Semi automatic operation
- Complete with pouch, battery + 1 set of electrodes
- Lightweight : 490g - Fits into a pocket

R 20 750.00 Excl VAT

Available units : 5

Price valid while stock last



SCHILLER

The Art of Diagnostics



EASY PULSE

- Latest Generation of Automatic CPR
- Efficient and compact mechanical CPR
- Unbelievably small and light electrically driven CPR device, weighing 3.5 kg
- Ideal for limited spaces
- Comprehensive
- Multidirectional chest compressions (3D)
- High-quality chest compressions of consistent rate and depth
- Possibility of 30:2 compression-ventilation cycle
- Easy to position – easy to operate
- Additional features
- Replaceable battery with charge level indicator
- Autonomous operation for 45 minutes (with one battery)
- Connections: external DC input, USB

Less
25%



R 158 850.00 Excl VAT

Available units : 2

Price valid while stock last



Physiogard Touch 7 Monitor



- Monitor for transport use
- ECG (12 lead), SpO2, NIBP
- 7" high definition touch screen device
- Touch 'n Save user interface
- 6 leads ECG monitoring
- Temperature (sensor to be ordered separately)
- Working on a Li-ion battery

R 132 000.00 Excl VAT including
Ambulance bracket SLIM
with DC/DC charging
module

Available units : 2

Price valid while stock last



Defigard Touch 7 Monitor/Defib



- Monitor / Defibrillator for transport use
- ECG (12 lead), SpO2, NIBP
- 7" high definition touch screen device
- Touch'n Save user interface
- AED defibrillator up to 200 J
- 6 leads ECG monitoring
- Temperature (sensor to be ordered separately)
- Working on a Li-ion battery

R 186 500.00 Excl VAT including
Ambulance bracket SLIM
with DC/DC charging
module

Available units : 1

Price valid while stock last



Defigard HD-7 Defibrillator



R 125 950.00 Excl VAT

- DEFIGARD HD-7 defibrillator / monitor with integrated printer, wifi and touch screen allowing :
- Manual and AED Defibrillation through the biphasic Multipulse Biowave® technology up to 200J
- With Paddles, AED, Pacer, 4-Lead ECG, Nellcor SpO2 and NIBP
- Large 7" Touch Screen with Touch 'n Save LCD interface
- A one-step mode selector to easily select the desired mode: semi-automatic AED, manual defibrillation, monitoring or pacing
- Integrated Thermal Printer

Available units : 1

Price valid while stock last



SCHILLER

The Art of Diagnostics



Defigard 400 Defibrillator



- DEFIGARD 400 defibrillator / monitor with integrated printer
- Large 8.4" colour TFT LCD screen with LED backlight
- Rotary knob for convenient input & on screen navigation
- Direct function keys for quick access
- Visual instructions and guidance on screen
- Stores Defibrillation Events
- 24 hours of Graphical & tabular trend for all parameters
- Adult and paediatric application
- Lightweight & portable

R 58 500.00 Excl VAT

Available units : 1

Price valid while stock last



SCHILLER
The Art of Diagnostics



DIAGNOSTIC STATION DS20

Less
25%



R 114 500.00 Excl VAT

- SCHILLER's DIAGNOSTIC STATION DS20
- Main vital signs and physical assessment tools in one device
- Large interactive touch screen 18.5" and self-explanatory user interface for maximum user-friendliness
- 12-Lead Resting ECG with Interpretation
- Masimo Pulse Oximetry,
- NIBP,
- Temperature

Available units : 1

Price valid while stock last



SCHILLER

The Art of Diagnostics



TRUSCOPE Vital Signs Monitor

Less
25%



- 7 inch TFT LCD Screen
- NIBP, SpO2, PR
- Rechargeable Lithium battery
- SpO2:
- Display of oxygen saturation and pulse rate
- Smart tones for pulse beep
- NIBP:
- Measurement principle: Oscillometric
- Auto, Manual and STAT Mode

R 37 450.00 Excl VAT

Including Mobile trolley

Available units : 13

Price valid while stock last

STAR8000 Patient Monitor



R 13 585.00 Excl VAT Unit Only

R 15 950.00 Excl VAT

Including Wall mount bracket

R 18 100.00 Excl VAT

Including Mobile trolley

- Multi Parameter patient monitor
- 12.1" TFT LCD Screen
- Standard: 5-lead ECG
- Comen SpO2
- NIBP
- TEMPERATURE
- RESP
- PR

Available units : 1

Price valid while stock last

C50 Patient Monitor



- Multi Parameter patient monitor
- 12.1" TFT LCD Screen
- Standard: 5-lead ECG
- Comen SpO2
- NIBP
- TEMPERATURE
- RESP
- PR

R 22 313.00 Excl VAT Unit Only

R 24 313.00 Excl VAT

Including Wall mount bracket

Available units : 14

Price valid while stock last

C80 Patient Monitor



- Multi Parameter patient monitor
- 15" TFT LCD Screen
- Standard: 5-lead ECG
- Comen SpO2
- NIBP
- TEMPERATURE
- RESP
- PR

Less
17%

R 30 200.00 Excl VAT Unit Only

R 32 550.00 Excl VAT

Including Wall mount bracket

Less
19%

Available units : 1

Price valid while stock last

S8 Defibrillator



Less
18%

R 72 520.00 Excl VAT

- 8.4 inch screen, provide 25 types energy selections from 1 ~ 360 J;
- Excellent battery performance, which supports 420 times 360 J release (with optional 2 batteries installed);
- 80mm width standard printer;
- Greater range of impedance range (15~250Ω) enables the best performance on various patients;
- High capacity and voice recording facilitate users to review the whole process of rescue.
- Operating Mode: Defibrillation, Monitoring, Pacing and AED mode
- Standard: 3/5-Lead ECG, Defibrillation Paddles, One Battery

Available units : 1

Price valid while stock last